# CFX Manager

#### Run set up

Protocol must be defined BEFORE run please check <u>your kit manual</u> for guidelines

Plate Layout can be edited after the run

# Run set up

	Bio-Rad Cl	X Manage	3.1 (ad	min)					
F	File View	User Run	Tools	Windows	Help				
	24 NAT	3		<b>X</b>					BIO RAD
	投 Detected li	nstrument(s)							
						Startup Wizard Run setup Repeat run Analyze	Select instrument Select run type	CFX96	
1. select	t us	er-	de	fine	ed ru	un			
		rument(c)							User:admin 10/16/2014 11:30 AM

# Run set up: protocol

2. make your own protocol or select	existing	3. you can edit existing protocol
File       View       User       Run         Image: Comparison of the state o	L (edmin )         Tools W ndows Help         Image: Create New         Select Existing         Selected Protocol         CFX_2stepAmp.prcl         Preview         Est. Run Time: 01:09:00 (96 Wells-All Channel         95.0 C         3:00	Express Load Express Load Ex
No instrument detected after 120	seconds.	User:admin 10/16/2014 11:32 AM 🤮

#### Run set up: protocol



#### Run set up: protocol



7. make your own plate layout o existing	or se	elect	t				8. y	/ou (	can (	edit	exist	ing I	ayou	ıt
existing	I (admi Tools V Run Set Run Set C Selectr QuickF Preview Fluorop A B C D E F G H	D) Vindows FUP otocol Create New ed Plate Plate_96 wells Vohores: 1 Unk Unk Unk Unk Unk Unk Unk	Help Plate All Channe FAM, H 2 Unk Unk Unk Unk Unk Unk	Start Run Start	ed, Cy5, Quas 4 Unk Unk Unk Unk Unk Unk Unk	sar 705 5 Unk Unk Unk Unk Unk Unk Unk	6 Unk Unk Unk Unk Unk Unk	7 Unk Unk Unk Unk Unk Unk	Plate Type 8 Unk Unk Unk Unk Unk Unk	Express Load QuickPlate_9 BR Clear 9 Unk Unk Unk Unk Unk Unk Unk	10         0         10         Unk         Unk	nannels.pltd Edit Se Scan Mode 11 Unk Unk Unk Unk Unk	ELEVENCE BIOR ELEVENCE All Channels 12 Unk Unk Unk Unk Unk Unk Unk Unk Unk	
No instrument detected after 120	seconds.										Use	er:admin 10	)/16/2014 11:3	8 AM

Eil	je Dav	Editor -	New	(admin)								9.	type	in your targets and templates in the experiment settings
	File	Settings	Editing T	ools 5can Mode	All Channels	~	😤 wel	l Groups	Trace	Styles	💷 Spreads	sheet View/I	mporter	Plate Loading Guide
	A B C D E F G		2	3	4	5	6	7	8	9				Sample Type   Load   Target Name   SYBR   Coad   Sample Name   Coad   Coad   Sample Name   Coad   Replicate #   1   Replicate Series
No ir	Pla	te Type: BR	Clear Vi	ew Sample	Well Group	b 🗌 Biolo	gical Set	] Well Note	3					Experiment Settings.         Clear Replicate #         Clear Wells         OK         OK         Cancel         User:admin       10/16/2014 11:40 AM

	12. and do the same for your templates	10. you can remove target names from the list
Bis Dad CEV Hanson 2 1   File   Plate Editor - New   File   Settings   Editing   100%   1   2   A   1   2   A   B   C   D   E   F   G   H	Image: Select To Arrow Select T	reet View/Importer
	New:     Add     Re       ○ Show Analysis Settings     Exclude the following sample types from Gene Expression     11. Or       ✓ NTC     NRT     Negative Control     Posit	<pre>checked item(s)</pre> Clear Wells Clear Wells Clear Wells OK Cancel
Plate Type: BR Clear	View Sample 🔲 Well Group 🔲 Biological Set 📄 Well Note econds.	OK Cancel User:admin 10/16/2014 11:45 AM:

	12. and do the same for your templates	10. you can remove target names from the list
Bis Dad CEV Hanson 2 1   File   File   Settings   Editor   100%   1   2   2   1   2   2   2   3   3   4   1   2   4   1   2   4   1   2   4   1   2   4   1   2   4   1   2   4   1   2   4   1   2   2   4	Image: Select To Arrow Select T	reet View/Importer
	New:     Add     Re       ○ Show Analysis Settings     Exclude the following sample types from Gene Expression     11. Or       ✓ NTC     NRT     Negative Control     Posit	<pre>checked item(s)</pre> Clear Wells Clear Wells Clear Wells OK Cancel
Plate Type: BR Clear	View Sample 🔲 Well Group 🔲 Biological Set 📄 Well Note econds.	OK Cancel User:admin 10/16/2014 11:45 AM:



14. select the wells and assign a sample type for them Settings Editing Tools File D H 😫 Plate Loading Guide 🔯 Scan Mode 🛛 All Channels Trace Styles... 💷 Spreadsheet View/Importer 100%  $\mathbf{v}$ Well Groups... 69 1 2 3 4 5 6 7 8 10 11 12 Select Fluorophores... А Sample Type В Unknown andarc Load NTC С Positive Control Negative Control NRT D Sample Name Load Е <none> 1 F Load Replicate # G Replicate Series Н -Experiment Settings... ź Clear Replicate # **e** Clear Wells View 0K Plate Type: BR Clear Cancel Sample 🔄 Well Group 📃 Biological Set 📃 Well Note No instrument detected after 120 seconds. User:admin | 10/16/2014 11:52 AM

15. select the wells and assign the target gene for them (you can type the name here directly)

File Settings Editing Tools														
	100%		5can Mode	All Channels		🕡 We	II Groups.	Trace	Styles	➡ Spreadst	heet View/Ir	nporter		Plate Loading Guid
A	1 Std	2 Std	3 Std	4	5	6	7	8	-	10	11	12	Select	Fluorophores
в	Std	Std	Std										Sample Type	Standard
с	Std Std	Std Std	Std Std										Load	Target Name
D E	Std	Std	Std										Load Sa	Gene 2 GOI Reference 1 Reference 2
F	Std	Std	Std										Load Re	plicate #
G													E 1	licate Series
н													Experir	nent Settings
														r Replicate #
Plate	Type: BR		ew Torget Mar			View	<b>—</b>						ОК	Cancel



18. in dilution series you can type in initial concentration (for primer efficiency use a random number) and the dilutions factor. and click on Apply



20. in the settings you can define what units you are going to use for your standards

	1 Plai	te Size	•	R/FAM	only	👌 We	I Groups	💐 Trace	Styles	💷 Spreadst	ieet View/Im	porter			Plate Loading	g Guide
	Nur	mber Conver	ntion	4	5	6	7	8	9	10	11	12		Select El	uorophores	
А	Uni sample A	ts L sample A	L sample A	<b>~</b> 0	opy number											
	Std-2	Std-2	Std-2	1 	old dilution	- 8							Sample 1	Гуре		~
В	1.00E+05 sample A	1.00E+05 sample A	1.00E+05 sample A	г	anomoles								Load		Target Name	9
С	<b>Std-3</b> 1.00E+04 sample A	Std-3 1.00E+04 sample A	<b>Std-3</b> 1.00E+04 sample A	F f	picomoles Temtomoles								SYB	R	<none></none>	~
D	<b>Std-4</b> 1.00E+03 sample A	<b>Std-4</b> 1.00E+03 sample A	<b>Std-4</b> 1.00E+03 sample A	r	attomoles nilligrams								Load	Sam Knor	iple Name ne>	~
E	<b>Std-5</b> 1.00E +02 sample A	<b>Std-5</b> 1.00E +02 sample A	<b>Std-5</b> 1.00E +02 sample A	r	nicrograms nanograms	-1							Load	Biol	ogical Set Nam	ie
F	<b>Std-6</b> 1.00E+01 sample A	<b>Std-6</b> 1.00E+01 sample A	<b>Std-6</b> 1.00E +01 sample A	F f	emtograms									<no< td=""><td>ne&gt;</td><td>~</td></no<>	ne>	~
G				F	attograms percent									Hep 1	licate #	× v
														Replic	ate Series	
Н														Experime	ent Settings	
													ß	Clear F	Replicate #	
													<b>1</b>	Cle	ar Wells	

21. open the lid, place your plate on the thermoblock , close the lid and start run



please make sure, that you have spun down your plate before running the reaction.

1. open your file (you can also drag and drop it into the software window) File Run Tools Windows Help User  $\sim$ BIO RAD 24 ? 🗙 Open Data File Detected Instrument(s) X G 🕸 📂 🖽-Look in: 🞯 Desktop My Documents G 🔜 My Computer 🧐 My Network Places My Recent Documents 🔗 VMware Shared Folders 📝 Data 2014-10-15 1434 Anna B Type: Bio-Rad iQ5 Data File Desktop Date Modified: 10/15/2014 4:09 PM Size: 919 KB Ô My Documents ID/Bar Code olume My Computer File name: Open My Network Files of type: All Data Files (\*.pcrd,\*.mgxd,\*.gxd\*,\*.opd,\*.tad,\* 🔽 Cancel Select All Blocks 19 1 Flash Block Indicator Open Lid Close Lid ≽ Start Run << Prev Next>> No instrument detected after 120 seconds. User:admin 10/16/2014 3:55 PM













go to the quantification data to see more information, you can sort the whole table be clicking on the small arrows in the header row. For example you can sort all data according to the target name



on the graph you will see the standard curve, Efficiency, R<sup>2</sup>, slope and y intercept values

please remind yourself, what are the desired and what are the allowed values for Efficiency and R<sup>2</sup>





#### Here you run Relative qPCR analysis



📶 Data .	Analysis - ad	lmin_2014-10-15	14-31-	28_BR002605	Zaenab.pcr	d									
File V	iew Settings	Export Tools	- <mark></mark> -	1elt Curve 🛃 1	Melt Curve Data	Gene	Expression	👤 End Point 🤞	Custom Data 1	view 🔁	QC 🔊	IP • 😱 Target Run Information	7	1 Jack	
111	Relative Normalized Expression			Stability Value	e Target 1 2 3		nt Variance 0.8921 1.0580 0.7956	M Value	3.8736 2.4663 2.2036			Mode: Normalized expression Graph Data: Relative to control Control Sample: None Show Chart Settings Show Chart Settings Experiment Se Target Stability	n (ΔΔCq) Value	A CONTRACTOR OF A CONTRACTOR O	1. Click on the Target stability value
	Target 4	△ Sample	♦ Cu	Recommended st Homogeneous Heterogeneous	Mean CV ability values fo CV < 0.25 < 0.5	0.9152 r samples: M < 0.5 < 1	Mean DEM	M Value 2.8475	ок	_	×		Ā		please remind yourself, what
	Reference 1 Reference 1 Reference 2 Reference 2 Reference 2 Reference 3 Reference 3	DNA pool sample 1 sample 2 DNA pool sample 1 sample 2 DNA pool sample 1			N/A N/A N/A N/A N/A N/A N/A N/A	N/A N/A N/A N/A N/A N/A N/A	N/A N/A N/A N/A N/A N/A N/A	32.22 36.79 30.81 31.31 27.51 28.25 36.78 24.40	0.88622 1.63328 0.41733 0.71607 0.18159 0.14935 0.72794 0.2542				=	And And And And	is it for and what are the desirable values
Completed	Reference 3	sample 1 sample 2 sample 2	can Mode:	All Channels Pla	N/A	N/A N/A ear Analysis	N/A N/A Mode: Baseline S	34.40 35.05 Subtracted Curve	0.23042 0.16716 Fit			User:admin 10/16	5/2014 10:23 PM		

2. adjust your selections in the Experiment Settings based on the M values you get